


The Contour handpiece by ENDYMED

# Redefining body treatments

The new Contour handpiece by EndyMed combines the power of focused **3DEEP energy** with the benefits of tissue lifting to provide the ultimate in cellulite reduction and body contouring treatments. The Contour's tissue lifting effect doubles the energy delivery and volumetric heating of the treatment area, achieving higher treatment temperatures for superior contouring, tightening and cellulite reduction. Improved blood circulation allows for better cell metabolism and the elimination of fat.

## Contour handpiece features:

- Unique electronically controlled modulated vacuum (ECMV™) doubles the deep heating volume for increased treatment effectiveness
- New modulated massage function for optimal cellulite reduction treatments
- Precise power and vacuum digital controls are located on the handpiece itself, for easy treatment experience
- Real time temperature tissue impedance monitoring - the best in the industry
- Completely safe and painless treatments, comfortable for the patient and operator



Contour deep thermal image: heating depth of more than 30mm

**ENDYMED**™  
3DEEP SKIN SCIENCE

\* The Contour treatment is FDA-cleared for temporary reduction in the appearance of cellulite and temporary improvement of local blood circulation

